


BOROSILICATE DISTILLATION GLASS (DOUBLE & SINGLE)



KEY FEATURES

- Made of borosilicate glass for high thermal and chemical resistance.
- Transparent setup allows clear visibility of distillation and condensation.
- Cooling Water from the condenser is re circulated into the boiling chamber hence saving energy
- Equipped with standard ground glass joints for leakproof connections.
- The stand is powder coated for rust free operation and lustrous appearance.
- Drain nozzle/ Teflon Stop Cork is provided for easy cleaning and maintenance.

TECHNICAL MATRIX

Model	PTE DG-302 A	PTE DG-302 B	PTE DG-302 C	PTE DG-302 D
Output Capacity	2	2	4	4
Stage level	Single	Double	Single	Double
Load	1.5 KW	3.0 KW	3.0KW	6.0 KW
Material	Borosilicate Glass			
Specific Conductivity	0.8-1us/cm			
Distillate Quality	Pyrogen Free Water			
Organic Matter	NIL			
Supply	220/230 Volts Single Phase			