

DRY BATH INCUBATOR



Model	PTE-DB-323
MOC	MS Duly powder coated
Display	LCD Display
Temp Range	R.T.+ 5°C ~ 100
Timing Range	60 min
Temp. Control Accuracy	+1.0°C
Temperature Resolution	0.1°C
Temp. Uniformity	+2°C
Heating Time	≤ 30min (25°C up to 100°C)
Standard Block	24 x 1.0ml/2.0 ml
Power	300W
Power supply	AC 220V,50/60Hz

INFRARED LOOP STERILIZER



Model	PTE IS 324 A	PTE IS 324 B
Max-Diameter Of Steriliser	Ø14	Ø35
Length Of Heating Zone	150mm	100mm
Power	150W	280W
Max Temp.	825°C ±50°C	
Heating Temp Mode	Two 400°C & 850°C	
Heating Mode	2	
Warm-up-Time	10-20 min	
Fuse	250V 3A 5°20	
Ambient Temp	5°C to 40°C	
Power Supply	AC220V/50Hz	
Gross Weight	1kg	1kg