

## HOT PLATE (ROUND/RECTANGULAR)



### KEY FEATURES

#### Sturdy Construction:

- Made of powder-coated mild steel or stainless steel for long-lasting durability and corrosion resistance.

#### Top Heating Surface:

- Top plate made from cast iron or SS 304 ensures uniform heat distribution and efficient heating.

#### Uniform Temperature:

- Provides consistent and even temperature across the surface for reliable results.

#### Temperature Control:

- Equipped with energy regulator or digital PID temperature controller for accurate temperature settings.

#### Temperature Range:

- It Operates from ambient to 350°C / 400°C, suitable for various laboratory applications.

#### Efficient Heating Element:

- Uses nichrome wire and ceramic insulated heating elements for quick and uniform heating performance.

### AVAILABLE SIZES

Model No.	B x L x H	Load
PTE HP 305 A	250 X 300 X 150 mm	1000 Watt
PTE HP 305 B	250 X 400 X 150 mm	1200 Watt
PTE HP 305 C	300 X 450 X 150 mm	2000 Watt
PTE HP 305 D	450 X 600 X 150 mm	3000 Watt
PTE HP 305 E	8" round	1000 Watt

- **Note :** Temp range in Round hot Plate will be upto 250 °C Granite , Ceramic Top is also Available in Some Models