

MAGNETIC STIRRER

TECHNICAL MATRIX

Model	PTE-MS-306	PTE-MS-307
Stirring Capacity	2 litres	5 Litre
Speed Range	100 – 2000 RPM	
Hot Plate	Available with or Without	
Top Plate Material	Stainless Steel	
Temperature Range	Ambient above +5°C to 110°C	
Speed Control	By Knob	
Magnet Type	PTFE-coated magnetic bar	
Top Plate Material	Stainless Steel	
Motor Type	DC Motor	
Power Supply	230V AC, 50Hz	

MULTI POSITION MAGNETIC STIRRER



OPTIONAL - ADD ONS

- Digital Speed Indicator and Digital Timer
- Digital Temperature Display
- External Temperature Probe (Pt100)
- Magnetic Stirrer Bars
- Temp Upto 300 /400 C
- Top Ceramic & Imported Model
- Capacity upto 20 Litres available

TECHNICAL MATRIX

Model	PTE-MPMS-308
Construction	MS Duly powder coated
Capacity	Up to 2 Litres Only
Speed Range	200-2000 RPM
Hot Plate	Available with or Without
Top Plate	Stainless Steel
Temperature Range	Ambient above +5°C to 110 °C
Speed Control	By Knob
Magnet Type	PTFE-coated magnetic bar
Motor	DC Motor
Supply	220/230 Volts AC supply
Optional	Digital RPM Indicator and Digital Temp Controller