

ROTA MANTLE



KEY FEATURES

- 1. Uniform Heating:** It is designed to Provide heating while allowing rotation /stirring in the flask
- 2. Temperature Control:** Equipped with an energy regulator or digital temperature controller for precise heat management upto 300 °C
- 3. Durable Construction:** Outer body made of powder-coated mild steel for durability and heat resistance.
- 4. Compatibility :** It is used for Round Bottom flask upto capacity of 500 ml to 5 Litre
- 5. Replaceable Mantle:** Heating element knitted into glass yarn can be easily replaced when needed. (Nichrome Wire Heating Element)

Model	PTE -RM313 A	PTE-RM313 B	PTE-RM313 C	PTE-RM313 D	PTE-RM313 E
Capacity	500ml	1000 ml	2000 ml	3000 ml	5000 ml

KEY FEATURES

- **Uniform Heating :** A Heating Mantle is a static electrical heater used to uniformly heat round-bottom flasks without rotation or stirring.
- **Temperature Control:** Equipped with an energy regulator or digital temperature controller for precise heat management upto 300 °C
- **Durable Construction:** Outer body made of powder coated mild steel for durability and heat resistance.
- **Replacable Mantle:** Heating element knitted into glass yarn can be easily.

HEATING MANTLE



Model	PTE -HM314 A	PTE-HM314 B	PTE-HM314 C	PTE-HM314 D	PTE-HM314 E
Capacity	500ml	1000 ml	2000 ml	3000 ml	5000 ml