

VORTEX SHAKER/MIXTURE PTE-119



TECHNICAL SPECIFICATIONS:

Speed Range:	0 to 3,000 rpm
Orbital (mixing motion):	4.0 mm to 6.0 mm
Operating modes	Touch, Continuous
Vessel capacity	50 mL tube or Ø30 mm
Body material	Die-cast aluminium base,
Power consumption	12 W / ~15 W
Dimensions	130 × 185 × 145 mm
Weight	3.3 kg
Protection / IP rating	IP21, IP42
Electrical supply	220-240 V AC, 50/60 Hz

VDRL ROTATOR/ SHAKER PTE-120



TECHNICAL SPECIFICATIONS:

Platform size	300 × 300 mm (12" × 12")
Speed (Rotation/RPM):	70-210 rpm Variable
Orbit (mixing motion):	10 mm orbit
Timer function	: 0-60 minutes (analog/digital)
Dimensions (W×D×H) :	300 × 300 × 125 mm
Platform holder	: Spring-wire or holder clamp to hold slides/plates/flasks
Platform MOC	: Stable Platform Anti Skid
MOC	: Powder-coated sheet metal, With heavy base motor
Digital Timer	Available (Optional)
Digital RPM Indicator	Available (Optional)