

WRIST ACTION SHAKER



Model No	Capacity
PTE-WAS-117 A	To Hold 4 Flask
PTE-WAS-117 B	To Hold 8 Flask
PTE-WAS-117 C	To Hold 12 Flask
PTE-WAS-117 D	Top Tray with 6 clamps

KEY FEATURES

- **Unique Wrist-Action Motion:** Simulates the gentle shaking action of the human wrist, ensuring thorough and uniform mixing.
- **Adjustable Speed Control:** Variable speed settings (typically 50–300 oscillations per minute) for precise control over agitation intensity.
- **Variable Stroke Length:** Adjustable shaking angle or stroke (usually 0°–15° or more) to suit different sample types.
- **Multiple Flask Capacity:** Can hold several flasks, bottles, or test tubes simultaneously with adjustable clamps upto 50 ml to 1000 ml
- **Sturdy Construction:** Heavy-duty steel body with powder-coated finish for durability and corrosion resistance. (stainless steel Arms and Brass Holders)
- **Quiet and Vibration-Free Operation:** Balanced drive mechanism ensures smooth performance and minimal noise.
- **Continuous or Timed Operation:** Integrated digital or mechanical timer for preset shaking durations.
- **Safety Features:** Overload protection and non-slip platform or clamps for secure sample holding.

OPTIONAL ADD ONS

- DIGITAL TIMER
- DIGITAL RPM INDICATOR
- DC DRIVE
- STAINLESS STEEL BODY
- OTHER CUSTOMIZATIONS RELATED TO CAPACITY AS WELL